



sartorius

Mouse IgG Type and Titer Kit



intellicyt[®]

A SARTORIUS BRAND

Open immediately upon
arrival and store reagents
at temperatures stated
on labels.

For Research Use Only.

Notice to Purchaser

The Mouse IgG Type and Titer Kit is a member of the IntelliCyt product line that has been tested extensively for live cell analysis applications. These screening kits are validated as complete screening assays and are optimized for use in high content screening applications. IntelliCyt's building blocks and reagents are designed for flexibility in multiplexing and incorporation into screening assays. IntelliCyt reagent kits are specifically formatted for optimal performance on IntelliCyt Screening platforms.

This product is manufactured and sold by INTELLICYT CORPORATION for research use only. The kit and components are not intended for diagnostic or therapeutic use. Purchase of the product does not include any right or license to use, develop, or otherwise exploit this product commercially. Any commercial use, development or exploitation of this product without the express written authorization of INTELLICYT CORPORATION is strictly prohibited.

Limited Warranty

These products are offered under a limited warranty. The products are guaranteed to meet appropriate specifications described in the product insert at the time of shipment. IntelliCyt Corporation will provide product replacement for valid claims. All claims should be made within five (5) days of receipt of order.

Trademarks and Patents

"Screening Solutions for Life" is a trademark of IntelliCyt. IntelliCyt, HyperView, HTFC, iDM, iQue, MultiCyt, ForeCyt, QBeads and Cy-Clone are registered trademarks of IntelliCyt. IntelliCyt's product portfolio is covered by US and foreign issued and pending patents.

List of Catalog Numbers





Description	Catalog No
Mouse IgG Type and Titer Kit 1x384 wells	91165
Mouse IgG Type and Titer Kit 5x384 wells	91166
Mouse IgG Type and Titer Kit 1x96 wells	91168







Kit Contents

Component	Quantity Provided
Mouse IgG Capture Beads	1 vial
Mouse FITC-IgG	1 vial
Mouse IgG Standard	1 vial
FL4 Membrane Integrity Reagent	1 vial
Sample Reaction Buffer	1 bottle
BSA (lyophilized)	1 bottle

NOTE: Add Sample Reaction Buffer to BSA **before** use.

Detection Channels

iQue Screener Standard Detector Channels					
Detector	Spectrum	Blue Laser (488 nm)		Red Laser (640 nm)	
533/30 nm		FL1	Bead Detection Channel		
585/40 nm		FL2			
670 nm LP		FL3			
675/25 nm				FL4	Cell Viability Detection Channel

iQue Screener PLUS Detector Channels (VBR)								
Detector	Spectrum	Violet Laser (405 nm)		Blue Laser (488 nm)		Red Laser (640 nm)		
445/45 nm		VL1						
530/30 nm		VL2		BL1	Bead Detection Channel			
572/28 nm		VL3		BL2				
615/24 nm		VL4		BL3				
675/30 nm		VL5		BL4		RL1	Cell Viability Detection Channel	
780/60 nm		VL6		BL5		RL2		

NOTE: For iQue Screener PLUS with VYB lasers, the bead detection channel is BL1 and the cell viability detection channel is VL3. For iQue Screener PLUS with BR lasers, the bead detection channel and the cell viability detection channel are the same as iQue Screener PLUS with VBR lasers.

Materials Needed but Not Provided

- IntelliCyt™ iQue® Screener system
- ForeCyt® Software
- Centrifuge capable of spinning microcentrifuge tubes and/or 15 mL conical tubes at up to 8,000 g
- Centrifuge capable of spinning microplates
- Vortex mixer
- Fresh complete cell culture media (Same media used to grow your sample cell culture)
- Microcentrifuge tubes and/or 15 mL conical tubes
- Two 50mL reagent reservoirs (Example source: VWR, Cat#89094-680) for reagent preparation purposes
- Universal black lid (Example source: Corning, Cat#3935) or foil to protect from light/evaporation
- 12-channel pipette reservoir (Example source: VWR, Cat#80092-466), optional for preparing serial dilutions
- Appropriate liquid handler or multi-channel pipette (Example Source: BioHit, Cat#LH745441, Picus NxT, 12channel, 5 to 120 µl; See [Appendix E.](#))


Table of Contents

List of Catalog Numbers	i
Kit Contents	i
Detection Channels	ii
Materials Needed but Not Provided	iii
Table of Contents	iv
Introduction	1
Assay Principles	2
Mix-and-Read Assay Workflow Overview	3
Best Practices and Tips	4
Assay Protocol Overview	5
Before You Begin	7
Reagent Preparation	7
Assay Set Up	14
Sample Acquisition and Data Analysis	17
Visualization of Screening Results	23
Appendix A: High Sensitivity Assay Protocol for Mouse Primary B-cell Samples	24
Appendix B: FAQ	25
Appendix C: Mixing Samples with the IntelliCyt Shaker	28
Appendix D: Plate-type Recommendations and Automated Wash Protocols for Microplates	29
Appendix E: Liquid Handler Recommendations	30
Abbreviated List of Consumables for IntelliCyt Screening Systems	31

Introduction


The Mouse IgG Type and Titer Kit was designed for ease of use in multiplexing. It provides a straightforward no-wash, no-sample dilution workflow that doesn't compromise assay performance and offers these unique advantages:

- **Precision** – Multiplexed isotype-capture bead assay enables more precise quantitation of IgG for each isotype.
- **Simultaneous isotyping measurement** – Provides clone purity information when two or more isotypes are present in the same well, and determines the DNA primer set for downstream gene cloning.
- **Simultaneous measurement of cell count and cell viability** – Monitors cell proliferation and cell health in the original cell cloning plates.
- **No-sample dilution assay** – Wide dynamic range (0.05 µg/mL to 50 µg/mL) enables transfer of cell culture samples directly into assay plates without a dilution step.
- **No-wash assay** – Mix and read format minimizes screen time, cost, and variability.



Assay Principles


The Mouse IgG Type and Titer Kit is a no-wash assay that enables the simultaneous measurement of mouse IgG isotypes, the mouse IgG quantity for each isotype, cell number and cell viability from each well of the screening plates. This is a competition assay in which fluorescently-labeled mouse IgG (mouse FITC-IgG) is added to mouse samples containing secreted IgG from hybridomas or B-cell cultures. The mouse FITC-IgG and non-labeled mouse sample IgG compete for binding to IgG capture beads in an isotype-specific manner. The amount of IgG isotype present in the sample is inversely proportional to the isotype-specific bead-associated fluorescence. Signals across four different isotype-specific beads determine the IgG isotype in the assay well. Cell viability is measured simultaneously in each well using cell membrane integrity dyes— fluorescent molecules that are cell impermeant. Healthy cells with intact cell membranes exclude the dye and are not fluorescent. Unhealthy cells with compromised membranes allow entry of the dye into the intracellular space where it then localizes in the nucleus and binds to DNA by intercalation.



Mouse IgG Type and Titer Kit assay principles. The no-wash competition assay functions on the differential binding of cell-secreted IgG vs mouse FITC-IgG to IgG Capture Beads in an isotype-specific manner. Samples with high IgG concentrations will exhibit a low fluorescence signal, whereas the absence or low concentration of solution IgG would result in a large fluorescence signal. The assay dynamic range is between 0.05 $\mu\text{g/mL}$ to 50 $\mu\text{g/mL}$.

Quantitative readouts from this assay can be measured as fluorescence intensity or extrapolated to a concentration ($\mu\text{g/mL}$) in solution via the use of isotype-specific standard curve.

Mix-and-Read Assay Workflow Overview



Screening Workflow. Distribute mouse hybridoma fusion cells (grown for the appropriate time) into culture plates by either limiting dilution or single cell cloning or single cell sorting. After growth, transfer 20 μ L samples from each well (supernatant only or supernatant plus cells) to assay plates. Mix first with mouse FITC-IgG and FL4 cell membrane integrity dye in the assay plates, and then mix with mouse IgG capture beads. After incubation at room temperature for 60 minutes, read plates directly on the IntelliCyt iQue Screening System.

Best Practices and Tips

If You're Running the Assay in 96-well Format


The assay protocol described in this manual is designed for 384-well plate format. However, if you're performing the assay in 96-well plates, you may use the same protocol and volumes designed for a 384-well format but adjusted to a 96-well format according to the volume table we provide in this manual (Pages 8-9). IntelliCyt recommends the use of 96-well bottom plates (IntelliCyt, Cat#90151). This assay kit provides ForeCyt templates for both 384-well and 96-well formats.

Manual Pipetting Recommendation

This protocol requires pipetting 5 μL volumes of liquid. If you are pipetting manually instead of using an automatic liquid handler, be extra careful during the 5 μL volume transfer of the prepared reagent from the reservoir to the assay well. If the plate is empty, touch the tip to the well bottom and then release all the liquid to transfer 5 μL volume into the well. If the plate already has a reagent or sample in the wells, touch the pipette tip to the upper inner wall of the well at 45-degree angle before you release the 5 μL prepared reagent. Touching the wall of the well prevents the 5 μL liquid droplet from hanging on the pipette tip instead of releasing into the assay well. A five second spin in a plate centrifuge will force the prepared detection reagent to the well bottom to mix with the existing reagent/sample already in the well.

Shaking (for Quick Mixing)

This assay requires shaking the micro-titer plate to quickly mix the sample/reagents. If you don't have a separate shaker, you can use the one on the IntelliCyt iQue screening system. (1) Click on **Device** in the menu bar. (2) Scroll down to **Manual Control**. (3) In the Manual Control window, use the arrows to set the RPM to 2000. (4) As soon as you click **On**, the shaker will begin to shake and continue to shake until you unclick.



Steps for using the shaker on the iQue screener.

Assay Protocol Overview

This protocol is designed to measure mouse IgG quantity and isotype in **mouse hybridoma samples** containing 1–50 µg/mL IgG. If you plan to measure mouse IgG and isotype in **mouse primary B** cell samples containing less than 1 µg/mL IgG (0.1–2 µg/mL IgG), refer to [Appendix A](#) for instructions for a High sensitivity assay. The following diagram is an overview of the protocol to help you plan your work. Detailed protocol instructions are provided in the next section.

Reagent Preparation

1

Add Sample Reaction Buffer to BSA

2

Dilute Capture Beads [1:5]

3

Mix 2 Detection Dyes:
1:125 FITC IgG
1:200 Membrane Integrity

4

Prepare IgG Standard with
FRESH CULTURE MEDIA



Assay Set Up

1

Add **5 μ L/well** Detection Dye Mixture

2

Add **20 μ L/well** standards/samples
Spin (500 rpm, 5 seconds)
Shake (2000 rpm, 20 Seconds)

3

Add **5 μ L/well** diluted capture beads
Spin (500 rpm, 5 seconds)
Shake (2000 rpm, 20 seconds)

4

Incubate 1 hour at room temperature protected from light



Sample Acquisition Data Analysis

Before You Begin

- Briefly centrifuge all vials before use to prevent reagent loss.
- Mix the dye with pipette or briefly vortex prior to use.
- Vortex capture beads prior to use to ensure a homogenous solution throughout the procedure. Beads tend to settle and aggregate over time. As with all things in life, mix them up occasionally.

Reagent Preparation

1. Mix Sample Reaction Buffer with BSA

- Add the entire volume of **Sample Reaction Buffer** solution to the lyophilized bovine serum albumin (BSA) bottle. Label the BSA bottle "**Sample Reaction Buffer with BSA**".
- Gently mix by inverting the bottle multiple times until the BSA has completely dissolved. Let the bottle sit at room temperature until you are ready to use it. If there are still particles, gently mix again before use. BSA does not affect mouse IgG quantification.



2. Dilute Mouse IgG Capture Beads (1st Reservoir)

Label a 50 mL reservoir "**Capture Beads**" and follow these steps for 1 x 384 well assay:

- Add 1.92 mL of the Sample Reaction Buffer with BSA you prepared earlier to the reservoir.
- Vortex the Mouse IgG Capture Beads vial from the kit for 30 seconds.
- Transfer 480 μ L of the Mouse IgG Capture Beads to the Sample Reaction Buffer in the Capture Beads reservoir.
- Mix the beads in the buffer by manual pipetting or by other gentle agitation. The capture beads are now ready to use. However, be patient and **DO NOT** transfer the prepared capture beads to the assay plate just yet. Cover reservoir with foil to prevent evaporation and protect from light.



For other kit sizes, refer to the following table to prepare the reagents.


If you are running a 96-well plate format, a 384-well kit can be used to run four 96-well plates or you may directly use 1 x 96-well kit. The dilution formula for a 96-well plate is listed in the table so you don't have to do the math.

	Sample Reaction Buffer with BSA	Mouse IgG Capture Beads (1:5 dilution)
1 x 384 wells	1.92 mL	480 µL
5 x 384 wells	9.6 mL	2.4 mL
1 x 96 wells	0.48 mL	120 µL

The volumes defined above will create enough prepared dye to add 5 µL per well for a full plate with minimal overage. To prepare stain for partial plates, or, for more overage, dilute the reagents at the dilution factors specified in sample reaction buffer to the desired total volume.

3. Mix 2 Detection Dyes (2nd Reservoir)

Label the second 50 mL reservoir "**Detection Reagent**". For 1 x 384-well assay you will add two dyes, **Mouse FITC-IgG** and **FL4 Membrane Integrity** dye.



- To the **Detection Reagent** reservoir, add 2.4 mL of the Sample Reaction Buffer with BSA.
- Add 19.2 µL of **Mouse FITC-IgG** to the Sample Reaction Buffer with BSA in the Detection Reagent reservoir.
- Add 12 µL of **FL4 Membrane Integrity** dye to the Detection Reagent reservoir.
- Mix the dye reagents in the buffer by manual pipetting. The reagent is now ready for use. However, be patient and do **NOT** transfer the mixed dye to the assay plate just yet. Cover reservoir with foil to prevent evaporation and protect from light.

For other kit sizes, please refer to the following table to prepare the reagents. If you are running a 96-well plate format, a 1 x 384-well kit can run 4 96-well plates, or you may directly use 1 x 96-well kit. The dilution formula for a 1 x 96 wells plate is also provided in the table.

	Sample Reaction Buffer with BSA	Mouse FITC-IgG (1:125 dilution)	FL4 Membrane Integrity Dye Reagent (1:200 dilution)
1 x 384 wells	2.4 mL	19.2 µL	12 µL
5 x 384 wells	12 mL	96 µL	60 µL
1 x 96 wells	0.6 mL	4.8 µL	3 µL

The volumes defined above will create enough prepared dye to add 5 µL per well for a full plate with minimal overage. To prepare stain for partial plates or to increase overage, calculate the total volume required and dilute the Mouse FITC IgG 1:125 and FL4 Membrane Integrity Dye 1:200.

4. Prepare Mouse IgG Standard

The kit contains one vial of Mouse IgG Standard with 4 mouse IgG isotype proteins (a mixture of IgG1, 2a, 2b and 3 at a concentration of 200 µg/mL for each isotype). The Mouse IgG Standard is used to generate a standard curve between 0.05 µg/mL to 50 µg /mL for each isotype.


The following dilution protocol is a guideline you can adapt or optimize for your laboratory (see Figure on Page 11).

- For a 384-well plate, prepare 16 micro-tubes. Label them #1 – 16. (To dilute the standard you may also use an empty 96-well plate or 12-channel pipette reservoir.)
- To tubes # 2-16, add 50 µL of FRESH cell culture media. This is the same media used to grow your sample culture. After you have added the cell culture media, set aside.
- To tube #1, add 75 µL of FRESH cell culture media. Then add 25 µL of the stock mouse IgG standard from the kit to the media in tube #1. The mixture should total 100 µL. Mix well by pipetting up and down six times. This sample is the highest concentration sample (50 µg/mL/isotype) for the standard curve (top standard).


- ❑ From tube #1, remove 50 μL of standard, and transfer to tube #2. Gently pipet up and down at least 6 times to completely mix the solutions.
- ❑ From tube #2, transfer 50 μL volume to tube #3 and mix by pipetting. Continue transferring and mixing until you reach tube #15. Do not transfer any standard into tube #16. Tube 16 will be the media only negative control.

The stock mouse IgG standard contains 200 $\mu\text{g}/\text{mL}$ for each mouse IgG isotype. Following these instructions will yield a standard curve top concentration of 50 $\mu\text{g}/\text{mL}/\text{isotype}$. If your sample has IgG concentration higher than 50 $\mu\text{g}/\text{mL}$, you may NOT need an intermediate sample dilution step. Please refer to Appendix B (FAQ, Question 4) for more information. If you need more mouse IgG standard than provided in the kit you may obtain additional standard from IntelliCyt (Cat#91171, Cat#91172, Cat#91173). Please note: ForeCyt software has capability to apply the standard curves from one plate with standards to any plates without standards in data analysis.

For 384-well Plate: 16-point Serial Dilution of IgG Standard.



For 96-well Plate 12-point Serial Dilution of IgG Standard




DO NOT add the Mouse IgG standards into the plate yet. Close the tubes containing the Mouse IgG standard to prevent evaporation; they will be transferred to the assay wells later.

Plate Design for IgG Standards

Below is an example of how to fill wells with the Mouse IgG Standards.

This design is already included in the template provided in the kit. We recommend using 1-2 rows of standards and arranging the standard wells from top to bottom and from low concentration to high concentration in 384 well plate; and, from left to right and from low concentration to high concentration in a 96 well plate. In order to achieve this configuration, choose **Edit Standard Set** under **Standards** in the **Design** tab and choose **Reverse Series**.




Arrange the standard wells from left to right (from low concentration to high concentration) in 96-well plate and from top to bottom (from low concentration to high concentration) in 384-well plate.

The screenshot shows the 'Add Standard Set' dialog box. The 'Name' field contains 'Standard Set 1'. The 'Wells' field is set to 'A01 to H02 by Column'. The 'Reverse Series' checkbox is checked. The 'Color' field shows an orange circle. Under 'Dilution Settings', the 'Dilution Factor' is set to '2'. The 'Set lowest concentration to 0' checkbox is checked. An 'Apply' button is visible at the bottom right.

Reverse Series to achieve a left to right (from low concentration to high concentration) in 96-well plate and top to bottom and (from low concentration to high concentration) in 384-well plate. Check the Set lowest concentration to zero checkbox.

Don't fill wells with the IgG standards until you read the Assay Set Up.

Warning: When you add the standards into the wells as described in Assay Set Up, follow the standard sequence you specified in Design when you initially set up your experiment. The wrong sequence will make the IgG quantitative number uninterpretable!



Representative standard curves (mouse IgG1, 2a, 2b, and 3) with 1:2 serial dilution and the top concentration at 50 μ g/mL/isotype.

Assay Set Up

This assay uses a no-wash workflow and provides results in terms of IgG concentration (e.g., $\mu\text{g}/\text{mL}$). You will use the serially diluted **Mouse IgG Standards** to generate the 4 standard curves that are used to measure the concentration for each of the mouse IgG isotypes.

The Capture Bead vial in the kit has a mixture of 4 types of capture beads. Each capture bead type has a specific affinity for a single mouse IgG isotype. The FITC detection signal on the capture beads has an inverse relationship with the mouse IgG concentration. Total mouse IgG concentration will be calculated in each well by ForeCyt software by adding up the 4 mouse IgG isotype concentrations. Cell counts and cell viability will be also analyzed for each well, if the sample includes cells.

This assay protocol is designed to screen and measure a mouse hybridoma culture sample with an IgG concentration between 1 $\mu\text{g}/\text{mL}$ and 50 $\mu\text{g}/\text{mL}$. If your sample IgG range is between 1 $\mu\text{g}/\text{mL}$ and 200 $\mu\text{g}/\text{mL}$, refer to Appendix B (FAQ, Q4) for more information on how to run a no-wash, no sample dilution assay. If your samples are from a primary mouse B cell culture with an IgG concentration between 0.1 $\mu\text{g}/\text{mL}$ and 2 $\mu\text{g}/\text{mL}$, please refer to Appendix A for high sensitivity assay protocol.

Total Protocol Time: 60 minutes


Total Hands-On Time: Approximately 15 minutes

NOTE: The sequential steps in the assay set up instructions below are critical to the success of your experiment. After you have prepared the reagents, the assay components must be added in this exact order:

1. Add Detection Reagent Mixture (Mouse FITC-IgG and FL4 Membrane Integrity Dye), 5 $\mu\text{L}/\text{well}$.
2. Add IgG Sample/IgG Standard, 20 $\mu\text{L}/\text{well}$; Quick spin. Mix.
3. Add prepared Capture Beads, 5 $\mu\text{L}/\text{well}$; Quick spin. Mix. Leave covered and protected from light for 60 minutes while you get coffee or something.

1. Add Detection Reagent Mixture (Mouse FITC-IgG and FL4 Membrane Integrity Dye)

To the bottom of each well, add 5 μ L of the Detection Reagent mixture (mouse FITC-IgG and FL4 membrane integrity dye) from the Detection Reagent reservoir.



2. Add Mouse IgG Sample/IgG Standard

- Add 20 μ L of mouse IgG sample (either IgG supernatant, or IgG supernatant/cell mixture) to each well you designated as a Sample.

The 20 μ L volume is appropriate for mouse hybridoma samples with IgG concentration between 1–50 μ g/mL. If your samples contain between 1–200 μ g/mL of IgG, please refer to FAQ, Question 4, in Appendix B for more information about designing a no-wash, no-sample-dilution measurement.

- Transfer 20 μ L Mouse IgG Standard prepared earlier to each well designated for IgG Standards.
- After you finish adding mouse IgG samples and standards to the plate, do a quick spin of the plate (500 g, 5 seconds) to ensure that all solutions are at the well bottom and not attached to the well sides.
- Mix the plate using the plate shaker on the IntelliCyt screening system for 20 seconds at 2,000 RPM. This ensures thorough mixing.



3. Add the Prepared Capture Beads Mixture

- Briefly mix the pre-diluted capture beads in the Capture Beads reservoir by gentle pipetting to keep the beads in suspension. Add 5 μL of the pre-diluted capture beads mixture to each well. Mix the beads in the reservoir 1–2 times while preparing a full assay plate to prevent beads from precipitating to the reservoir/tube bottom.
- Do a quick spin of the plate (500 g, 5 seconds) to ensure that samples are at the bottom of the wells and not attached to the sides.
- Mix the plate to ensure thorough mixing. Use the shaker on the IntelliCyt screening system and mix for 20 seconds at 2,000 RPM.




4. Incubate at room temperature for 60 minutes with lid on the plate, protected from light

After incubation, the plate is ready for sample acquisition on the iQue screener. Make sure there is a QSol® buffer cartridge in the first rinse station to prime the tubing.




Sample Acquisition and Data Analysis

1. Launch ForeCyt screening software.
2. If you have not already imported the template for this assay on the USB drive in the kit, insert the USB drive in your computer. The following screenshot shows how to import the template.




3. Find and select the corresponding assay template for your iQue platform, software version and plate type (96w or 384w) to run this assay.
4. Create a New Experiment using the template. The following screenshot shows the steps.




5. The templated acquisition settings (found in the USB drive provided in the kit) are in bold and are found in the bulleted list below.
- Pre-Plate Prime: **60 seconds**. Prior to sampling a plate, the system will prime the tubing with the QSol buffer.
 - Pre-Plate Shake: **15 seconds at 2,400 RPM**. Prior to sampling a plate, the shaker will agitate the samples at the specified speed and time. The specified shake speed assumes that the sample contains 30 μ L of total sample volume. For samples that contain different volumes, see Appendix C.
 - Sampling Order: **By row**. The order of acquisition and how the probe moves from well to well can be specified by column or a zig-zag pattern.
 - Sip Time: **2 seconds**. Sip time is proportional to the sample volume collected per well. Sip time can be increased to acquire more events per well. The tradeoff is slightly slower overall read-time.
 - Additional Up Time: **0.5 seconds**. Up Time specifies the amount of time the probe pauses before moving to the next well and determines the spacing between samples. We recommend no adjustment.
 - Pump Speed: **Standard speed (15 RPM on iQue Standard, 29 RPM on iQue PLUS)**. This specifies the rate at which samples pass by the detectors. The default setting of the standard roughly corresponds to introduction of \sim 1.5 μ L per second. We recommend no adjustment. This approximate 1.5 μ L volume may not be precise. For absolute event counting, please refer to Appendix B (FAQ, Q11).

- **Interwell Shake: 4 seconds at 2,400 RPM after every 12 wells if sampling by row or after every 16 wells if sampling by column.** If you're using a 96-well plate, specify a shake speed of 2000 RPM. Inter-well Shake inserts a rinse and shake after a specified number of wells to keep sample particles in suspension. Adjust the frequency of Inter-well Shake if you find sample particles precipitating to the well bottom.
 - **Speed (for iQue Screener or HTFC): Medium.** Specify the flow rate of the detector. This affects the resolution of the data. We recommend no adjustment.
 - **Threshold: FSC-H at 80,000 on iQue Standard; 100,000 on iQue PLUS.** Specify the lower thresholds for data acquisition. The current setting on FSC will filter out events with size below 2 microns. We recommend no adjustment.
6. Once you've completed any Protocol adjustments, click **Run** on the **Controller** to acquire data.




7. During the plate read, data will automatically begin populating into the pre-defined analysis template.
8. Verify that the sample data aligns with the pre-defined gating strategy from the template. If necessary, adjust the gates in each plot to encompass the proper bead populations. All gates can be moved by clicking the gate label and dragging to the desired location.

Identify Cell and Bead Populations. If necessary, move the "cells" gate and "beads" gate to encompass all regions of interest. While the ForeCyt template will gate the cells and beads for you, it might miss the smaller bead population so you'll have to adjust the gate. The smaller bead population **must** be included in the gate of all beads. You can also enlarge the size of these gates if you want to include additional cell populations in the analysis.




Gate Beads and Cells.

Identify Singlet Beads. Analyzing only the singlet bead population helps avoid analysis artifacts created when analyzing aggregates of beads. The singlet population is on the $\sim 45^\circ$ angle on the FSC-H vs FSC-A plot.




Singlet bead gate.

Identify 4 Different Isotype-specific Capture Beads. All 4 mouse IgG isotype-specific capture beads from singlet beads must be gated as shown below. For iQue Screener PLUS, create a RL1-H vs FSC-H plot from singlet beads, and gate 4 isotype-specific populations. As with the other pre-gated plots, make any necessary adjustments.



Gate 4 Isotype-specific Capture Beads on iQue Screener PLUS


For iQue Screener, the gating strategy is slightly different as shown below. Create a FSC-A vs SSC-A plot with the singlet beads population. Gate mouse IgG2b capture beads first (the population on the right), and then gate other 3 beads. Next, create FL4-A vs FL3-A plot from "Other 3 beads" population, and gate mouse IgG1, IgG2a, IgG3 population from the left to the right.




Gate 4 Isotype-specific Capture Beads on iQue Screener.

Identify Viable Cells. Histograms show the measurement of cell viability in the RL1-H channel (iQue Screener PLUS) or the FL4-H channel (iQue Screener standard). This is a function of the membrane integrity dye binding to damaged cells' DNA. In the figure below, positive cells (those not encompassed in the Live Cell gate) represent the non-viable population. Adjust the gate (if needed) to encompass only the left-most peak. This represents the viable cells that excluded the dye. This gate separates the percentage of viable cells in each well.

NOTE: The FL4 membrane integrity dye used in the assay identifies both necrotic cells (highly stained population furthest to the right) and apoptotic cells (intermediate stained population between viable and necrotic cells).








Viable cells on RL1-H (iQue Screener Plus) gated in ForeCyt.



Viable cells FL4-H (iQue Screener Standard) gated in ForeCyt.

Visualization of Screening Results

After all the gates have been verified and adjusted as necessary for the plate-level data set, you can generate additional analyses including heat maps, dose responses, standard curves, and Panoramas.




Data visualization examples.



Appendix A: High Sensitivity Assay Protocol for Mouse Primary B-cell Samples

Normal hybridoma cultures contain IgG concentrations ranging from 1 to 50 $\mu\text{g/mL}$. In contrast, mouse primary B cell cultures, the IgG concentration ranges from 0.1 to 2 $\mu\text{g/mL}$. The following protocol increases the assay sensitivity allowing accurate quantification of mouse IgG in B-cell cultures:

1. Add 5 μL prepared capture beads to each well
2. Add 20 μL prepared sample or standards to each well. Quick spin (500 g, 5 seconds). Briefly mix the plate (2,000 rpm, 20 seconds)
3. Incubate at room temperature for 120 minutes without light
4. Add 5 μL prepared 2 detection reagent mixture to each well. Quick spin (500 g, 5 seconds). Briefly mix the plate (2,000 rpm, 20 seconds)
5. Incubate at room temperature for 60 minutes without light. Now the plate is ready for sampling on iQue Screener platform.



Appendix B: FAQ

Q1: Why does the standard curve look like a straight flat line at saturation level instead of an S curve?

Answer: When there is a straight line instead of a standard curve at saturation level, the competition assay setup is wrong. You must mix the mouse FITC-IgG with mouse IgG standard first. Add the capture beads at the end. You must strictly follow the competition setup protocol to run this assay. The assay standard curve should look like a reverse S shape, not a normal S shape. When IgG concentration goes higher, the bead detection fluorescence goes lower. If your standard curve is a normal S curve, not a reverse S curve, it is very likely your standard design is opposite. Please go to ForeCyt design, click "standard" tab, and then click "Edit Standard Set", and finally check or uncheck the box "Reverse Series". If you need to change standard design in multiple plates, you may use "Copy Plate Design" function and make sure to check the box "standards" in the pop-up window.

Q2: In a multi-plate screening campaign, is there any fluorescence signal shifting on capture beads?

Answer: This assay is a no wash competition assay. After assay reaction, the fluorescence signal on the beads is stable for 2-3 hours. After 3 hours, there may be a slight increase of fluorescence signal compared with the signal at Time 0. If you have standards assigned in each plate, there will be no issue even after 6 hours because the IgG quantitation will be based on the standard curves created from the same plate and this will not change the IgG quantitation for unknown samples. However, if you only have one standard plate and do not assign standard wells in each screening plate, you may re-sample the standard plate after every 2 hours of acquisition, and then apply the standard curves generated at each specific time point to the subsequent screening plates within 2 hours acquisition time in ForeCyt software.

Q3: Does this assay work with IgG quantification for other species such as rat?

Answer: No. This assay will only measure IgG from mouse species; no other species will work with this assay.

Q4: Do I need to dilute my sample if my sample IgG range is between 1-200 µg/mL?

Answer: This is a no-wash no-dilution assay for measuring mouse IgG from hybridoma samples with most common range between 1-50 µg/mL. However, if you expect your sample IgG might sometimes be within 1-200 µg/mL range, there are 2 options to handle this issue in order to maintain the no-wash, no-sample-dilution workflow:

Option 1: You may use 5 µL volume for each sample instead of 20 µL volume at the assay step when you add the sample into the well. Then, add 15 µL fresh complete culture media to each sample well. In this add-mix workflow, you dilute your sample 4 fold but avoid the tedious intermediate dilution-and-then-transfer step. For the standard wells, you still use 20 µL volume. After IgG quantitation with standard curve, the output IgG quantity will be 4 times lower than its actual quantity, and you may use ForeCyt software advanced metrics function to multiply the calculated sample IgG quantity by 4 to get actual sample IgG quantity.

Option 2: The recommended top concentration for the standard is 50 µg/mL/isotype in the protocol. However, the kit provides the standard stock at 200 µg/mL/isotype. You may use 200 µg/mL for the standard as the top concentration. The volume of mouse IgG is limited, so be conservative when preparing the dilution series for the standard curves. It is very important to use just 5 µL volume for both unknown sample and standard. In this way, the total reaction volume is 15 µL: 5µL detection reagent + 5 µL sample or standard + 5 µL capture beads. Test to make sure your plate type can accommodate 15 µL sample acquisition. Under most circumstances there is no problem in sample acquisition with a 15 µL volume for 384 well plate. However, for 96 well plate, you may need to confirm that the geometry of your plate can handle the 15µL volume in the sample acquisition on IntelliCyt screening platform. If there is a problem, you can add an additional 15 µL sample reaction buffer with BSA to each well after the assay reaction is finished, and just before sample acquisition on the IntelliCyt screening system. You may double your sip time per well if necessary in the Protocol tab of ForeCyt.

Q5: I have bead number variation in the assay plate. Does bead number variation decrease my assay performance?

Answer: Bead variation generally does NOT affect assay performance at all, as long as you have acquired a hundred or more total beads per well. The bead number variation does NOT change the median fluorescence intensity (MFI) of the capture beads. The IgG quantification is based on bead MFI. You can avoid bead number variation by agitating the diluted beads in reservoir periodically during the manual transferring of prepared beads to assay plate. You may also increase the sip time in the Protocol tab of ForeCyt before you sample the whole plate.

Q6: I have difficulty in precisely pipetting 5 µL volume in the assay plate. How should I better handle the transferring of 5 µL volume from reservoir to assay plate?

Answer: An automated liquid handler can pipette 5 µL volume very precisely. If you do not use an automated liquid handler, we recommend a 12-channel pipette (5-120 µL) for the liquid transfer (see Appendix E). The pipette tip must touch the wall of the well at 45-degree angle before you manually release the 5 µL prepared reagent. Touching the inner wall of the well prevents the 5 µL liquid droplet from hanging on the pipette tip instead of releasing into the assay well. A quick spin will force the solutions into the bottom of the well.

Q7: My culture plate is a 96-well format. Can the 384-well assay kit be adapted to 96-well format?

Answer: Yes. A single 384-well assay kit can be used for 4 assay plates in a 96-well format or you may directly use 1 x 96-well kit to run 96-well plate assay. If performing the assay in 96-well plates, you may use the same protocol and protocol volumes designed for a 384-well format. IntelliCyt recommends the use of 96-well V-bottom plates (IntelliCyt, Cat#90151). You may need to adjust shake speeds and other protocol settings to accommodate the 96-well plate format. The assay kit provides assay templates for both 384-well format and 96-well format on the USB drive included in the kit package.

Q8: Should I clean my instrument before a screening campaign?

Answer: Routinely cleaning and maintaining the IntelliCyt screening system is very important for accurate results of your assay. Follow the hardware manual cleaning directions, regardless of the samples you run through the system. Make sure to use an IntelliCyt QSol cartridge on the buffer station (first rinse station) for tubing priming. In addition, do not use your own adhesive seal to seal the assay plate or the cartridge on any rinse station. Adhesive seals may block/clog the sampling probe or even the cytometer engine if the sampling probe touches your adhesive seal during the sampling process. If your cell samples do cause system clogging, change the sample tubing and fluidic link (if you are using iQue Screener PLUS). Do a long clean and an unclog by following software instruction in the ForeCyt Controller. Normally the clogging will go away. If not, repeat the long clean and unclogging processes.

Q9: Can I multiplex this assay with additional cell health endpoints or with antigen binding measurements?

Answer: Yes, it is possible. You may multiplex this assay with apoptosis measurement by using one or more dyes from IntelliCyt MultiCyt 4-plex Apoptosis Kit. You may also multiplex this assay with antigen binding measurement by using a fluorescent antigen or by using an unlabeled antigen plus a fluorescent antibody against your antigen. However, you may need to do extensive assay optimization of the workflow to confirm performance in your biology.

Q10: Can I fix my samples in the plate with PFA?

Answer: Yes, it is possible. You may add 10 μ L 4% PFA to each well after the assay reaction and mix the plate well by shaking at 2,000 rpm for 20 seconds. After 10 minutes incubation, you may sample the plate on the iQue Screener platform. The fixation may extend the detection signal on the capture beads. You should validate this against your biology and your plate type materials. Some plate type materials may cross-link with biology samples with fixatives. The cell viability information will not be precise since the fixation may break the cell membrane and allow the membrane dye to stain the nuclei of healthy cells.

Q11: Can I do absolute calculation of cell density in my sample?

Answer: It is possible to improve precision in cell density calculation. We recommend using SPHERO™ AccuCount Particles (Spherotech, Cat#ACBP-50-10) to do the volumetric measurement on the iQue platform. This bead has the absolute count per volume unit.

Please follow the Spherotech protocol to mix the beads well and transfer the beads to a testing plate. We recommend running the exact sample protocol and the same plate type and the same volume in the well as used in this Mouse IgG Type and Titer assay. You may only need to run 2-3 wells at the beginning of your assay day to measure the sip volume on your iQue platform. Only use this volume measurement on the same day of the experiment in which you wish to calculate the cell density. Please adjust your final calculation by considering the sip time (in the ForeCyt Protocol) and the dilution factor of your sample in the final assay reaction volume. In order to avoid cell count variance as a result of possible cell re-attachment to the well bottom/wall, we recommend using the Greiner cell-repellent plate (e.g., Cat#651970, Cat#781970) or Corning ultra-low binding plate (e.g., Cat#7007, Cat#4516) for a more precise cell count.

Note: The cell viability analysis based on fluorescence dye is uses a different method than traditional Trypan Blue-based viable cell measurement. The FL4 membrane integrity dye used in this assay not only stains the necrotic cells but also the apoptotic cells. The cell viability number may be lower in this assay than a Trypan Blue-based assay.

Appendix C: Mixing Samples with the IntelliCyt Shaker

iQue Screener and iQue Screener PLUS

Plate Type	Well Volume	MAX rpm
96-Well	20–40 μL	2,600
96-Well	40–60 μL	2,200
96-Well	60+ μL	A/O*
384-Well	10–30 μL	3,000
384-Well	30–50 μL	2,800
384-Well	50+ μL	A/O*

HTFC Screening System

Plate Type	Well Volume	MAX rpm
96-Well	20–40 μL	2,800
96-Well	40–60 μL	2,400
96-Well	60+ μL	A/O*
384-Well	10–30 μL	3,500
384-Well	30–50 μL	3,000
384-Well	50+ μL	A/O*

iQue Screener HD

Plate Type	Well Volume	MAX rpm
96-Well	20–40 μL	3,200
96-Well	40–60 μL	2,400
96-Well	60+ μL	A/O*
384-Well	10–30 μL	3,500
384-Well	30–50 μL	3,100
384-Well	50+ μL	A/O*
1536-Well	up to 6 μL	5,000

*A/O = Additional Optimization necessary. While it is possible to run these volumes, they were not routinely tested during assay development. To determine ideal shake speeds for high volume assays, IntelliCyt recommends starting at low RPM values and slowly increasing to higher values.

Appendix D: Plate-type Recommendations and Automated Wash Protocols for Microplates

The following plate types and aspiration settings have been extensively tested with IntelliCyt assay products.

Plate Type	Well Type	Manufacturer	Manufacturer Product #
384-well	V-bottom	Greiner	781280
96-well	V-bottom	IntelliCyt	90151

When using the above plate types, the following aspiration programs have been tested on a BioTek ELx405 Select. If you have a different plate washer brand or model, it is possible to approximate the aspiration settings on a different system.

IntelliCyt highly recommends that wash protocols use an automated plate washer. Manual aspiration of plates and/or plate inversion could result in severe sample loss.

Plate Type	Aspiration Height Setting	Aspiration Height Offset	Aspiration Rate Setting	Aspiration Rate
384-well, V-bottom	#31	3.937 mm	#6	15 mm/sec
96-well, V-bottom	#40	5.08 mm	#6	15 mm/sec

Appendix E: Liquid Handler Recommendations

IntelliCyt recommends the following liquid handlers:

12-channel pipette for manual transfer of liquid to the plate:

- Manual 12-channel pipette mLINE or Tacta, 5–100 μL (Sartorius);
- Electronic 12-channel pipette eLINE or Picus, 5–120 μL (Sartorius).

Automated liquid handler:

- Personal Pipettor, 96 or 384 channels (Apricot Designs);
- MicroFlow Select, 8 channels (BioTek).

Single-channel pipette for reagent preparation:

- Manual single-channel pipette mLINE or Tacta (Sartorius);
- Electronic single-channel pipette eLine or Picus (Sartorius).

Abbreviated List of Consumables for IntelliCyt Screening Systems

iQue®/iQue® Screener PLUS and HTFC® Probes

Part	Description
90659	iQue® Probe Et Tubing Assy for Gen 2 iQue® Et iQue® HD – 5 Pk
91088	iQue® Screener PLUS Probe Et Tubing Assy for iQue® Screener PLUS – 5 Pk
91093	iQue® Screener PLUS FluidLink tubing connector – 5 Pk

iQue®/iQue® Screener PLUS and HTFC® Solutions

Part	Description
90077	Decontamination Concentrate Solution for iQue® Screener/HTFC® – 5 PK (makes 1 Liter)
90078	Bacteriostatic Concentrate Soln for Sheath Fluid for iQue®/HTFC®/iQue® Screener PLUS
90079	Cleaning Concentrate Solution (makes 1 Liter)
90082	Extended Flow Cell Cleaning Solution for iQue®/HTFC®
90083	QSol Buffer Cartridge-Fluidic Station (Single) for iQue® Screener/PLUS/HD
90286	iQue® Fluidic Station Buffer Cartridge – 10 Pk
90287	QSol Buffer Cartridge-Fluidic Station (10 Pk) for iQue® Screener/PLUS/HD
90288	iQue® Fluidic Station Decon/Cleaner Cartridge – 10 Pk
90289	iQue® Fluidic Station Water Cartridge – 10 Pk
91089	PLUS ONE Detector Maintenance Solution only for iQue® Screener PLUS (orange label) 5X
91090	PLUS TWO Detector Maintenance Solution only for iQue® Screener PLUS (purple label) 5X
91304	QSol Buffer Concentrate Solution (100x; makes 500 mL) – Use in Sampling Area Fluid Station for iQue® Screener/PLUS/HD/HTFC® – Use directly in assay media/sample buffer for iQue® Screener/HD/HTFC®

iQue®/iQue® Screener PLUS and HTFC® Marker Beads

Part	Description
90040	FL1 In-Well Marker Beads for iQue®/HTFC® – 10 X 384 well plates
90041	FL2 In-Well Marker Beads for iQue®/HTFC® – 10 X 384 well plates
90042	FL3 In-Well Marker Beads for iQue®/HTFC® – 10 X 384 well plates
90043	FL4 In-Well Marker Beads for iQue®/HTFC® – 10 X 384-well plates
90044	In-Well Marker Beads for iQue®/HTFC® Starter Kit (4 colors)
90635	FL1 Between-Well Marker Cartridge for iQue® (ForeCyt 4.0 or later req'd)
90636	FL2 Between-Well Marker Cartridge for iQue® (ForeCyt 4.0 or later req'd)
90637	FL3 Between-Well Marker Cartridge for iQue® (ForeCyt 4.0 or later req'd)
90638	FL4 Between-Well Marker Cartridge for iQue® (ForeCyt 4.0 or later req'd)

iQue®/iQue® Screener PLUS and HTFC® Maintenance

Part	Description
90075	iQue®/HTFC® Fluidics Maintenance Kit
90295	6 peak Validation beads (for Red Laser and FL4 Detector)
90296	8 peak Validation Beads (for Blue Laser and FL1,FL2,FL3 Detectors)
91091	PLUS Validation Beads (all channels) for iQue Screener PLUS
91094	iQue® Screener PLUS Maintenance Kit
91095	iQue® Screener PLUS Maintenance Kit – 10 pack

**Refer to Intellicyt.com for complete list. Contact your local area sales representative for part number and pricing information

ABOUT THE COVER:

The New Mexico Thistle (*Cirsium neomexicanum*) grows at elevations up to 6500 feet throughout the Southwest. It can reach heights of up to 6 feet and blooms following winter and spring rains and again after heavy summer monsoon rainfall. It favors plains, mesas, rocky hillsides, foothills, roadsides and washes. New Mexico Thistle has been used for medicinal purposes by Navajo and Yavapai Indians and provides resources for bees, butterflies and other insects. Monarch butterflies visit New Mexico Thistle more than any other wildflower during their migration back to Mexico.

intellicyt[®]

A SARTORIUS BRAND

5700 Pasadena Ave. NE
Albuquerque, NM 87113
USA

T (505) 345-9075
F (866) 782-3140
E support.intellicyt@sartorius.com
www.intellicyt.com

Part Number: 11377, Rev A, March 2018