

Analyte Specifications:

Catalog Number:

Basic FGF


92012

MultiCyt™ QBeads PlexScreen


Reconstituting the protein Standard in 2 mL of assay diluent results in a concentration of 5000 pg/mL. Standard curves for no-wash and wash protocols were run in triplicate (+/- S.D.) starting from 5000 pg/mL and proceeding in a 1:2 serial titration. The limit of detection (dashed line) is 3 standard deviations above the mean of samples without Standard (N=12). Results are typical and may vary based on specific assay conditions.

STANDARD CURVES

No Wash Protocol



1-Wash Protocol



Lowest Detectable Test Conc. (pg/mL): 78.1

Lowest Detectable Test Conc. (pg/mL): 78.1

For research purposes only. Not for resale.