

Analyte Specifications:

Catalog Number:

**IL-21**



**92071**

**MultiCyt® QBeads™  
PlexScreen Mouse Secreted  
Protein**

Reconstituting the IL-21 standard in 2 mL of assay diluent results in a concentration of 20,000 pg/mL. Standard curves for no-wash and 1-wash protocols are shown starting from 20,000 pg/mL and proceeding in a 1:2 serial titration. The limit of detection, represented by the dashed line, is determined by adding 3 standard deviations to the mean of samples without IL-21 (N=12). Results are typical and may vary based on specific assay conditions.

**No-Wash Assay Performance**

**1-Wash Assay Performance**



**Lowest Detectable  
Concentration Tested (pg/mL): 10**

**Lowest Detectable  
Concentration Tested (pg/mL): 10**

For research purposes only. Not for resale.

Rev. 1 Sep 2014