

Analyte Specifications:

Catalog Number:

IFN γ



92082

**MultiCyt[®] QBeads[™]
PlexScreen Rat Secreted
Protein**

Reconstituting the IFN γ standard in 2 mL of assay diluent results in a concentration of 20,000 pg/mL. Standard curves for no-wash and 1-wash protocols are shown starting from 20,000 pg/mL and proceeding in a 1:2 serial titration. The limit of detection, represented by the dashed line, is determined by adding 3 standard deviations to the mean of samples without IFN γ (N=12). Results are typical and may vary based on specific assay conditions.

No-Wash Assay Performance

1-Wash Assay Performance



Lowest Detectable

Concentration Tested (pg/mL): 156

Lowest Detectable

Concentration Tested (pg/mL): 9.8

For research purposes only. Not for resale.

Rev. 1

Sep 2014